

# Running the gel

## 1. Prepare the samples

Make up your samples to 20  $\mu$ L by adding 10  $\mu$ L distilled water to each of the tubes using a micropipette.

Using a quick wrist flicking motion get all the liquid to the bottom of the tube so that you can pipette it out.

## 2. Prepare your E-gel-EX

Take the gel out of the foil packaging.

Carefully remove the comb by lifting it from both sides without bending it.

## 3. Load the gel

Pipette your samples and ladder into the desired wells; try to avoid introducing bubbles into the wells.

Fill any empty wells with 20  $\mu$ L of distilled water.

**Write down which well you loaded each sample.**

## 4. Run the gel

Place the gel on the grey iBase, right side first, sliding it across and pressing on the left side.

Use programme 7 "**Run E-Gel EX**" programme (find this by pressing Mode then the up arrow then press **Go** when you are ready).

## 5. Watch progress

You can watch the progress of the gel using the Safe-Imager and the orange screen: place the iBase on top of the Safe-Imager with the orange screen over the gel, then press the red button on the Safe-Imager to turn on the light – your DNA bands will light up.

## 6. View results and photograph

When the gels are finished (10 minutes), have a look at each of them using the Safe-Imager, and take photos if desired.

The DNA will diffuse over time, so be sure to examine your gels shortly after they finish running.

*If you have any problems call the invitrogen helpline:  
0800 335 997*



Add 10  $\mu$ L of distilled water to all 5 of the eppendorf tubes using a micropipette. Then transfer the solutions to the gel wells.



Safe Imager

iBase

Orange filter screen



Combined unit ready to run gel.